

PRAME RABBIT MAB

货号: N262751

产品全名: PRAME 兔单克隆抗体

基因符号 MAPE; OIP4; CT130; OIP-4

UNIPROT ID: P78395

背景: Functions as a transcriptional repressor, inhibiting the signaling of retinoic acid through the retinoic acid receptors RARA, RARB and RARG. Prevents retinoic acid-induced cell proliferation arrest, differentiation and apoptosis.

抗原: Recombinant protein of human PRAME

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-7B6

分子量: Calculated MW: 58 kDa; Observed MW: 58 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Mouse


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Tags & Cell Markers

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PRAME in K562 lysates using PRAME antibody.



Immunohistochemistry analysis of paraffin-embedded mouse testis using PRAME antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.