

PREB RABBIT MAB

货号: N261621

产品全名: PREB 兔单克隆抗体

基因符号 Prolactin regulatory element-binding protein; Mammalian guanine nucleotide exchange factor mSec12

UNIPROT ID: Q9HCU5

背景: Guanine nucleotide exchange factor that specifically activates the small GTPase SAR1B.

抗原: Recombinant protein of human PREB

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-8B7

分子量: Calculated MW: 45 kDa; Observed MW: 45 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PREB antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PREB in K562, C6, 3T3, HeLa lysates using PREB antibody.