

## PRESENILIN 2 RABBIT MAB

货号: N261566

产品全名: Presenilin 2 兔单克隆抗体

基因符号 AD3LP; AD5; E5-1; STM-2

**UNIPROT ID:** P49810

背景: Probable catalytic subunit of the gamma-secretase complex, an endoprotease complex that catalyzes the intramembrane cleavage of integral membrane proteins such as Notch receptors and APP (amyloid-beta precursor protein).

抗原: A synthetic peptide of human Presenilin 2

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-5B8

分子量: Calculated MW: 50 kDa; Observed MW: 23 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience&Alzheimer's disease

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Presenilin 2 in rat Brain, C6, 3T3, HeLa lysates using Presenilin 2 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---