

PRION PROTEIN RABBIT PAB

货号: N225667

产品全名: Prion Protein 兔多抗

基因符号 CJD; GSS; PrP; PRNP; ASCR; KURU; PRIP; PrPc; CD230; AltPrP;
Prion protein; PrP27-30; PrP33-35C

UNIPROT ID: P04156

背景: The PRNP gene encodes the major prion protein (PrP, CD230), a widely-expressed glycoprotein expressed at high levels in the central nervous system. While the typical cellular function of PrP is not well defined, it is a putative antioxidant and a metal-binding protein that may be involved in signal transduction.

抗原: A synthesized peptide derived from human Prion 蛋白

经过测试的应用: WB,IHC-P,ICC/IF,FC

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 28 kDa; Observed MW: 28 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience&Prion disease

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Prion Protein in Human fetal brain lysates using Prion Protein antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
