

PRPF4 RABBIT MAB

货号: N262771

产品全名: PRPF4 兔单克隆抗体

基因符号 PRP4; RP70; HPRP4; Prp4p; HPRP4P; SNRNP60

UNIPROT ID: O43172

背景: Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.

抗原: A synthetic peptide of human PRPF4

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-6C5

分子量: Calculated MW: 58 kDa; Observed MW: 58 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PRPF4 in K562, 3T3, HeLa lysates using PRPF4 antibody.