

PSD93 RABBIT MAB

货号: N263511

产品全名: PSD93 兔单克隆抗体

基因符号 DLG2; Disks large homolog 2; Channel-associated protein of synapse-110; Chapsyn-110; Postsynaptic density protein PSD-93

UNIPROT ID: Q15700

背景: Required for perception of chronic pain through NMDA receptor signaling. Regulates surface expression of NMDA receptors in dorsal horn neurons of the spinal cord. Interacts with the cytoplasmic tail of NMDA receptor subunits as well as inward rectifying potassium channels.

抗原: A synthetic peptide of human PSD93

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-4E7

分子量: Calculated MW: 98 kDa; Observed MW: 110 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PSD93 in HeLa, C6, Ramos lysates using PSD93 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
