

## PTF1A RABBIT MAB

货号: N261612

产品全名: PTF1A 兔单克隆抗体

基因符号 Pancreas-specific transcription factor 1a; bHLH transcription factor p48; PTF1-p48

**UNIPROT ID:** Q9QX98

背景: Transcription factor implicated in the cell fate determination in various organs. Binds to the E-box consensus sequence 5'-CANNTG-3'. Plays a role in early and late pancreas development and differentiation. Important for determining whether cells allocated to the pancreatic buds continue towards pancreatic organogenesis or revert back to duodenal fates.

抗原: Recombinant protein of mouse PTF1A

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-7B1

分子量: Calculated MW: 35 kDa; Observed MW: 48 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of PTF1A in K562, C6 lysates using PTF1A antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---