

R CADHERIN RABBIT MAB

货号: N261692

产品全名: R Cadherin 兔单克隆抗体

基因符号 cadherin 4; type 1; R-cadherin (retinal); CAD4; RCAD; R-CAD

UNIPROT ID: P55283

背景: Cadherins are calcium-dependent cell adhesion proteins.

抗原: A synthetic peptide of human R Cadherin

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-5H6

分子量: Calculated MW: 100 kDa; Observed MW: 130 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Tags & Cell Markers

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of R-Cadherin/Cadherin 4 in rat Brain, C6 lysates using R-Cadherin/Cadherin 4 antibody.



Western blot analysis of R-Cadherin/Cadherin 4 in mouse brain, U87-MG lysates using R-Cadherin antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
