

## RAB29 RABBIT MAB

货号: N263518

产品全名: Rab29 兔单克隆抗体

基因符号 RAB7L; RAB7L1

**UNIPROT ID:** O14966

背景: The small GTPases Rab are key regulators in vesicle trafficking.

抗原: Recombinant protein of human RAB29

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-6A3

分子量: Calculated MW: 23 kDa; Observed MW: 23 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of RAB29 in U251, U87-MG lysates using Rab29 antibody.