

## RANK RABBIT PAB

货号: N225136

产品全名: RANK 兔多抗

基因符号 CD265; FEO; LOH18CR1; mRANK; ODFR; OFE; OPTB7; OST5; PDB2; RANK; TNFRSF11A; TRANCER

**UNIPROT ID:** Q9Y6Q6

背景: Receptor for TNFSF11/RANKL/TRANCE/OPGL; essential for RANKL-mediated osteoclastogenesis. Involved in the regulation of interactions between T-cells and dendritic cells.

抗原: Synthetic peptide from human protein at AA range: 60-120

经过测试的应用: WB,IHC-P,ELISA

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 66 kDa; Observed MW: 66 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human liver using RANK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of RANK in Sodium Hela lysates using RANK antibody.