

RASSF1 RABBIT PAB

货号: N225861

产品全名: RASSF1 兔多抗

基因符号 RASSF1; RDA32; Ras association domain-containing protein 1

UNIPROT ID: Q9NS23

背景: Potential tumor suppressor. Required for death receptor-dependent apoptosis. Mediates activation of STK3/MST2 and STK4/MST1 during Fas-induced apoptosis by preventing their dephosphorylation. When associated with MOAPI, promotes BAX conformational change and translocation to mitochondrial membranes in response to TNF and TNFSF10 stimulation.

抗原: A synthesized peptide derived from human RASSF1

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 39 kDa; Observed MW: 39 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of RASSF1 in HeLa lysates using RASSF1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
