

## RB RABBIT MAB

货号: N261697

产品全名: Rb 兔单克隆抗体

基因符号 p110-RB1; pRb; pp105

**UNIPROT ID:** P13405

背景: Tumor suppressor that is a key regulator of the G1/S transition of the cell cycle.

抗原: Recombinant protein of mouse Rb

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-8G7

分子量: Calculated MW: 105 kDa; Observed MW: 105 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Rb in C6, Jurkat, HeLa lysates using Rb antibody.