

RNA HELICASE A RABBIT PAB

货号: N225964

产品全名: RNA Helicase A 兔多抗

基因符号 LKP; RHA; DDX9; NDH2; NDH11

UNIPROT ID: Q08211

背景: Unwinds double-stranded DNA and RNA in a 3' to 5' direction. Alteration of secondary structure may subsequently influence interactions with proteins or other nucleic acids. Functions as a transcriptional activator. Component of the CRD-mediated complex that promotes MYC mRNA stability.

抗原: A synthesized peptide derived from human RNA Helicase A

经过测试的应用: WB,IHC-P,ICC/IF,FC

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 141 kDa; Observed MW: 141 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of RNA Helicase A in Jurkat lysates using RNA Helicase A antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
