

RPA70 RABBIT MAB

货号: N263534

产品全名: RPA70 兔单克隆抗体

基因符号 RPA1; HSSB; MST075; REPA1; RF-A; RP-A; RPA70

UNIPROT ID: P27694

背景: Plays an essential role in several cellular processes in DNA metabolism including replication, recombination and DNA repair. Binds and subsequently stabilizes single-stranded DNA intermediates and thus prevents complementary DNA from reannealing.

抗原: A synthetic peptide of human RPA70

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-3F4

分子量: Calculated MW: 68 kDa; Observed MW: 68 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of RPA70 in A549, HL-60 lysates using RPA70 antibody.