

RRM2 RABBIT MAB

货号: N262851

产品全名: RRM2 兔单克隆抗体

基因符号 R2; RR2M

UNIPROT ID: P31350

背景: Provides the precursors necessary for DNA synthesis. Catalyzes the biosynthesis of deoxyribonucleotides from the corresponding ribonucleotides. Inhibits Wnt signaling.

抗原: Recombinant protein of human RRM2

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-4G8

分子量: Calculated MW: 45 kDa; Observed MW: 45 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of RRM2 in K562, HeLa lysates using RRM2 antibody.



Immunocytochemistry analysis of RRM2 (green) in Jurkat using RRM2 antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
