

RUNX RABBIT MAB

货号: N262857

产品全名: RUNX 兔单克隆抗体

基因符号 RUNX1; AML1; RUNX2; OSF2; RUNX3; CBFA3

UNIPROT ID: Q01196/Q13761/Q13950

背景: RUNX2 regulates the transcription of various genes including osteopontin, bone sialoprotein, and osteocalcin via binding to the core site of the enhancers or promoters. RUNX3/AML2 is a member of the Runt family of transcription factors. RUNX3 is important for the suppression of cell proliferation in the gastric epithelium, neurogenesis of the dorsal root ganglia, and T cell differentiation.

抗原: A synthetic peptide of human RUNX1/2/3

经过测试的应用: WB, IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-719

分子量: Calculated MW: 49 kDa; Observed MW: 49 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of RUNX1/2/3 in RAW264.7 lysates using RUNX1/2/3 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse colon using RUNX antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.