

## S100 ALPHA9 RABBIT MAB

货号: N262860

产品全名: S100 alpha9 兔单克隆抗体

基因符号 S100 alpha9; Calgranulin-B; Leukocyte L1 complex heavy chain

**UNIPROT ID:** P31725

背景: S100A9 is a calcium- and zinc-binding protein which plays a prominent role in the regulation of inflammatory processes and immune response.

抗原: Recombinant protein of mouse S100A9

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-3C3

分子量: Calculated MW: 13 kDa; Observed MW: 13 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of S100A9 in mouse spleen, mouse lung lysates using S100 alpha9 antibody.