

SALL4 RABBIT MAB

货号: N262866

产品全名: Sall4 兔单克隆抗体

基因符号 Zinc finger protein 797; Zinc finger protein SALL4

UNIPROT ID: Q9UJQ4

背景: Transcription factor with a key role in the maintenance and self-renewal of embryonic and hematopoietic stem cells.

抗原: Recombinant protein of human Sall4

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R02-6H7

分子量: Calculated MW: 112 kDa; Observed MW: 112 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Sall4 in Daudi lysates using Sall4 antibody.