

## SCAI RABBIT MAB

货号: N262870

产品全名: SCAI 兔单克隆抗体

基因符号 NET40; C9orf126

**UNIPROT ID:** Q8N9R8

背景: Tumor suppressor which functions to suppress MRTFA-induced SRF transcriptional activity. May function in the RHOA-DIAPH1 signal transduction pathway and regulate cell migration through transcriptional regulation of ITGB1.

抗原: A synthetic peptide of human SCAI

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-3J2

分子量: Calculated MW: 70 kDa; Observed MW: 60 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SCAI antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of SCAI in C6, 3T3 lysates using SCAI antibody.