

SCGB1A1 RABBIT PAB

货号: S216873

产品全名: SCGB1A1 兔多抗

基因符号 UGB; UPI; CC10; CC16; CCSP; UP-1; CCPBP

UNIPROT ID: P11684 (Gene Accession - BC004481)

背景: This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.

抗原: Fusion protein of human SCGB1A1

经过测试的应用: ELISA, IHC

推荐稀释比: IHC: 50-200; ELISA: 5000-10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

亚型: Immunogen-specific rabbit IgG


纯化: Antigen affinity purification

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using 216873(SCGB1A1 Antibody) at a dilution of 1/80(Cytoplasm).