

SLC18A3 RABBIT PAB

货号: S222847

产品全名: SLC18A3 兔多抗

基因符号 VACHT

UNIPROT ID: Q16572 (Gene Accession - NP_003046)

背景: This gene is a member of the vesicular amine transporter family. The encoded transmembrane protein transports acetylcholine into secretory vesicles for release into the extracellular space. Acetylcholine transport utilizes a proton gradient established by a vacuolar ATPase. This gene is located within the first intron of the choline acetyltransferase gene.

抗原: Synthetic peptide of human SLC18A3

经过测试的应用: ELISA, IF

推荐稀释比: ELISA: 5000-60000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

亚型: Immunogen-specific rabbit IgG

纯化: Antigen affinity purification

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing