

## SLC27A4 RABBIT MAB

货号: N262907

产品全名: SLC27A4 兔单克隆抗体

基因符号 SLC27A4; ACSVL4; FATP4; Long-chain fatty acid transport protein 4; FATP-4; Fatty acid transport protein 4; Solute carrier family 27 member 4

**UNIPROT ID:** Q6P1M0

背景: plays a role in the transport of long chain fatty acids across the plasma membrane. It has acyl-coA ligase activity for long chain and very long chain fatty acids.

抗原: Recombinant protein of human SLC27A4/FATP4

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-3C7

分子量: Calculated MW: 72 kDa; Observed MW: 72 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of SLC27A4/FATP4 (green) in HeLa using SLC27A4/FATP4 antibody, and DAPI (blue)



Western blot analysis of SLC27A4/FATP4 in K562, HeLa lysates using SLC27A4 antibody.