

SLIT2 RABBIT PAB

货号: N225873

产品全名: Slit2 兔多抗

基因符号 slit guidance ligand 2; SLIL3; Slit-2

UNIPROT ID: O94813

背景: Thought to act as molecular guidance cue in cellular migration, and function appears to be mediated by interaction with roundabout homolog receptors.

抗原: A synthesized peptide derived from human Slit2

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 170 kDa; Observed MW: 58 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Slit2 in SH-SY5Y lysates using Slit2 antibody.