

SMARCA5 RABBIT PAB

货号: N225984

产品全名: SMARCA5 兔多抗

基因符号 ISWI; SNF2H; HISWI; hSNF2H; WCRF135

UNIPROT ID: O60264

背景: Helicase that possesses intrinsic ATP-dependent nucleosome-remodeling activity. Complexes containing SMARCA5 are capable of forming ordered nucleosome arrays on chromatin; this may require intact histone H4 tails.

抗原: A synthesized peptide derived from human SNF2H

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 122 kDa; Observed MW: 122 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of SNF2H in MCF-7 lysates using SMARCA5 antibody.