

SNX4 RABBIT MAB

货号: N262922

产品全名: SNX4 兔单克隆抗体

基因符号 ATG24B

UNIPROT ID: O95219

背景: Involved in the regulation of endocytosis and in several stages of intracellular trafficking (PubMed:12668730). Plays a role in recycling endocytosed transferrin receptor and prevent its degradation.

抗原: A synthetic peptide of human SNX4

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-8A9

分子量: Calculated MW: 52 kDa; Observed MW: 52 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of SNX4 in C6, 3T3, Hela lysates using SNX4 antibody.