

SODD RABBIT MAB

货号: N262897

产品全名: SODD 兔单克隆抗体

基因符号 BAG4; SODD; BAG family molecular chaperone regulator 4; BAG-4; Bcl-2-associated athanogene 4; Silencer of death domains

UNIPROT ID: O95429

背景: Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release. Prevents constitutive TNFRSF1A signaling. Negative regulator of PRKN translocation to damaged mitochondria.

抗原: A synthetic peptide of human Silencer of Death Domain

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-6J7

分子量: Calculated MW: 50 kDa; Observed MW: 70 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Silencer of Death Domain in MOLT4 lysates using SODD antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
