

SOX17 RABBIT PAB

货号: N225160

产品全名: SOX17 兔多抗

基因符号 SOX17; Transcription factor SOX-17

UNIPROT ID: Q9H6I2

背景: Acts as transcription regulator that binds target promoter DNA and bends the DNA.

抗原: Synthesized peptide derived from Sox-17 . at AA range: 90-170

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性 : Rabbit Polyclonal

分子量: Calculated MW: 44 kDa; Observed MW: 44 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of SOX17 in Jurkat lysates using Sox17 antibody.