

SRP72 RABBIT MAB

货号: N262935

产品全名: SRP72 兔单克隆抗体

基因符号 BMFF; BMFS1; HEL103

UNIPROT ID: O76094

背景: Signal-recognition-particle assembly has a crucial role in targeting secretory proteins to the rough endoplasmic reticulum membrane. Binds the 7S RNA only in presence of SRP68. This ribonucleoprotein complex might interact directly with the docking protein in the ER membrane and possibly participate in the elongation arrest function.

抗原: A synthetic peptide of human SRP72

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R08-614

分子量: Calculated MW: 75 kDa; Observed MW: 75 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of SRP72 in HeLa, HUVEC, K562 lysates using SRP72 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
