

STAR RABBIT PAB

货号: N225884

产品全名: StAR 兔多抗

基因符号 STARD1

UNIPROT ID: P49675

背景: Steroidogenic acute regulatory protein (StAR) plays a significant role in cholesterol transport from the cytoplasmic outer membrane to the inner mitochondrial membrane.

抗原: A synthesized peptide derived from human StAR

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 32 kDa; Observed MW: 32 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of StAR in SKBR3 lysates using StAR antibody.