

STAT2 RABBIT PAB

货号: N226059

产品全名: STAT2 兔多抗

基因符号 Homo sapiens interferon alpha induced transcriptional activator; ISGF 3; P113; signal transducer and activator of transcription 2 113kD; STAT113; Stat2

UNIPROT ID: P52630

背景: STAT2 (113-kDa), originally purified from the nuclei of alpha-interferon-treated cells, is critical to the transcriptional responses induced by type I interferons, IFN-alpha/beta. Stat2 is rapidly activated by phosphorylation at Tyr690 in response to stimulation by IFN-alpha/beta via associations with receptor-bound Jak kinases

抗原: A synthetic peptide of human STAT2

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 98 kDa; Observed MW: 98 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of STAT2 in K562 lysates using STAT2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
