

STAT6 RABBIT MAB

货号: N263563

产品全名: STAT6 兔单克隆抗体

基因符号 STAT6; Signal transducer and activator of transcription 6; IL-4 Stat

UNIPROT ID: P42226

背景: STAT6 transcription factor of the STAT family. Plays a central role in IL4-mediated biological responses. Induces the expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4.

抗原: A synthetic peptide of human STAT6

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-0B5

分子量: Calculated MW: 94 kDa; Observed MW: 110 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of STAT6 in HeLa, U2OS lysates using STAT6 antibody.