

## STAT6 RABBIT PAB

货号: N226014

产品全名: STAT6 兔多抗

基因符号 STAT6; Signal transducer and activator of transcription 6; IL-4 Stat

**UNIPROT ID:** P42226

**背景:** STAT6 transcription factor of the STAT family. Plays a central role in IL4-mediated biological responses. Induces the expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4.

**抗原:** A synthetic peptide of human STAT6

**经过测试的应用:** WB,IHC-F,IHC-P,ICC/IF,FC,IP

**推荐稀释比:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 94 kDa; Observed MW: 110 kDa

**亚型:** IgG


**纯化:** Affinity Purified

**种属反应性:** Mouse

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Signal Transduction

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of STAT6 in Raw264.7 lysates using STAT6 antibody.