

## STING RABBIT PAB

货号: N225928

产品全名: STING 兔多抗

基因符号 ERIS; MITA; MPYS; SAVI; NET23; STING; hMITA; hSTING; STING-beta; Transmembrane Protein 173

**UNIPROT ID:** Q86WV6

**背景:** Facilitator of innate immune signaling that promotes the production of type I interferon (IFN-alpha and IFN-beta). Innate immune response is triggered in response to non-CpG double-stranded DNA from viruses and bacteria delivered to the cytoplasm.

**抗原:** A synthesized peptide derived from human TMEM173

**经过测试的应用:** WB, ICC/IF, FC

**推荐稀释比:** WB: 1/500-1/1000 IF: 1/50-1/200 FC: 1/50-1/100

**种属反应性:** Rabbit

**克隆性:** Rabbit Polyclonal

**分子量:** Calculated MW: 42 kDa; Observed MW: 37 kDa

**亚型:** IgG


**纯化:** Affinity Chromatography

**种属反应性:** Human

**成分:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**研究领域:** Immunology

**储存和运输:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TMEM173 in HeLa lysates using Transmembrane Protein 173 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---