

STIP1 RABBIT PAB

货号: N225924

产品全名: STIP1 兔多抗

基因符号 Hop; P60; STII; STIIL; STIP1; Stress induced phosphoprotein 1

UNIPROT ID: P31948

背景: Mediates the association of the molecular chaperones HSC70 and HSP90 (HSPCA and HSPCB).

抗原: A synthesized peptide derived from human STIP1

经过测试的应用: WB,IHC-P,ICC/IF,IP,FC

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 63 kDa; Observed MW: 63 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of STIP1 in Jurkat lysates using STIP1 antibody.