

STK33 RABBIT MAB

货号: N262947

产品全名: STK33 兔单克隆抗体

基因符号 STK33; Serine/threonine-protein kinase 33

UNIPROT ID: Q9BYT3

背景: Serine/threonine protein kinase which phosphorylates VIME. May play a specific role in the dynamic behavior of the intermediate filament cytoskeleton by phosphorylation of VIME (By similarity). Not essential for the survival of KRAS-dependent AML cell lines.

抗原: Recombinant protein of human STK33

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R03-9B3

分子量: Calculated MW: 58 kDa; Observed MW: 58 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of STK33 in Jurkat, 293 lysates using STK33 antibody.