

## STRIATIN 4 RABBIT MAB

货号: N262948

产品全名: Striatin 4 兔单克隆抗体

基因符号 ZIN; PPP2R6C; zinedin

**UNIPROT ID:** Q9NRL3

背景: Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein. Miscellaneous The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane).

抗原: Recombinant protein of human Striatin 4

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-4J1

分子量: Calculated MW: 81 kDa; Observed MW: 81 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Striatin 4 in HeLa lysates using Striatin 4 antibody.