

SURVIVIN (8B9) MOUSE MAB

货号: N261284

产品全名: Survivin (8B9) 小鼠单抗

基因符号 BIRC5; API4; IAP4; Baculoviral IAP repeat-containing protein 5; Apoptosis inhibitor 4; Apoptosis inhibitor survivin

UNIPROT ID: O15392

背景: Survivin is an apoptosis inhibitor that is expressed during the G2/M phase of the cell cycle. Associates with the microtubules of the mitotic spindle and any disruption results in the loss of apoptosis activity. May play a role in neoplasia.

抗原: Recombinant Protein of Survivin

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 8B9-6G6-1G4

分子量: Calculated MW: 16 kDa; Observed MW: 16 kDa

亚型: IgG1


纯化: Affinity Purified

种属反应性: Human,Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cell Biology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma using Survivin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Survivin in COS7, Jurkat and Hela lysates using Survivin antibody.



Western blot analysis of Survivin (8B9) in TF1, HL6 lysates using Survivin (8B9) antibody.