

SURVIVIN RABBIT MAB

货号: N262957

产品全名: Survivin 兔单克隆抗体

基因符号 BIRC5; API4; IAP4; Baculoviral IAP repeat-containing protein 5; Apoptosis inhibitor 4; Apoptosis inhibitor survivin

UNIPROT ID: O15392

背景: Survivin is an apoptosis inhibitor that is expressed during the G2/M phase of the cell cycle. Associates with the microtubules of the mitotic spindle and any disruption results in the loss of apoptosis activity. May play a role in neoplasia.

抗原: Recombinant protein of human Survivin

经过测试的应用: WB, ICC/IF, IP

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R04-1B1

分子量: Calculated MW: 16 kDa; Observed MW: 16 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics&Histone phosphorylation

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of Survivin (green) in MCF-7 using Survivin antibody, and DAPI (blue)



Western blot analysis of Survivin in Jurkat lysates using Survivin antibody.