

SWITCH ASSOCIATED PROTEIN 70 (2E9) MOUSE MAB

货号: N261075

产品全名: Switch Associated Protein 70 (2E9) 小鼠单抗

基因符号 HSPC321; SWAP-70; FLJ39540; KIAA0640

UNIPROT ID: Q9UH65

背景: Phosphatidylinositol 3,4,5-trisphosphate-dependent guanine nucleotide exchange factor (GEF) which, independently of RAS, transduces signals from tyrosine kinase receptors to RAC.

抗原: Purified recombinant human SWAP70 protein fragments expressed in E.coli.

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 2E9-E8-E9

分子量: Calculated MW: 69 kDa; Observed MW: 69 kDa

亚型: IgG2b


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat,Monkey


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunoprecipitation analysis of Switch Associated Protein 70 (2E9) in HeLa lysates using SWAP70 antibody.



Western blot analysis of SWAP70 in COS7, 293T, C6, 3T3 and HeLa lysates using SWAP70 antibody.