

T BOX PROTEIN 3 RABBIT MAB

货号: N262976

产品全名: T Box Protein 3 兔单克隆抗体

基因符号 T box 3; T-box protein 3; Tbx3; TBX3 ISO; UMS; XHL

UNIPROT ID: O15119

背景: Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.

抗原: A synthetic peptide of human Tbx3

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R07-5B9

分子量: Calculated MW: 79 kDa; Observed MW: 79 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Mouse

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Tbx3 in 3T3 lysates using T Box Protein 3 antibody.