

TACSTD2 RABBIT MAB

货号: N262964

产品全名: TACSTD2 兔单克隆抗体

基因符号 TACSTD2; GA733-1; MIS1; TROP2; Tumor-associated calcium signal transducer 2; Cell surface glycoprotein Trop-2; Membrane component chromosome 1 surface marker 1; Pancreatic carcinoma marker protein GA733-1

UNIPROT ID: P09758

背景: This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.

抗原: A synthetic peptide of human TACD2

经过测试的应用: WB,IHC-F,IHC-P,ICC/IF

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R05-312

分子量: Calculated MW: 36 kDa; Observed MW: 45-65 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cancer

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TROP2 in HCT116, PC-3, MCF-7 lysates using TACSTD2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
