

TAF3 RABBIT MAB

货号: N261572

产品全名: TAF3 兔单克隆抗体

基因符号 Transcription initiation factor TFIID subunit 3

UNIPROT ID: Q5VWG9

背景: The TFIID basal transcription factor complex plays a major role in the initiation of RNA polymerase II (Pol II)-dependent transcription.

抗原: Recombinant protein of human TAF3

经过测试的应用: WB,IHC-P

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-6A4

分子量: Calculated MW: 104 kDa; Observed MW: 140 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TAF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of TAF3 in Jurkat lysates using TAF3 antibody.