

TAGAP RABBIT MAB

货号: N262967

产品全名: TAGAP 兔单克隆抗体

基因符号 FKSG15; IDDM21; TAGAP1; ARHGAP47

UNIPROT ID: Q8N103

背景: May function as a GTPase-activating protein and may play important roles during T-cell activation.

抗原: Recombinant protein of human TAGAP

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R01-3D3

分子量: Calculated MW: 81 kDa; Observed MW: 81 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TAGAP in rat Brain, C6, CHO-K1, HeLa lysates using TAGAP antibody.