

TAOK1 RABBIT MAB

货号: N263570

产品全名: TAOK1 兔单克隆抗体

基因符号 PSK2; TAO1; KFC-B; MARKK; PSK-2; hKFC-B; hTAOK1; MAP3K16

UNIPROT ID: Q7L7X3

背景: Serine/threonine-protein kinase involved in various processes such as p38/MAPK14 stress-activated MAPK cascade, DNA damage response and regulation of cytoskeleton stability.

抗原: A synthetic peptide of human TAOK1

经过测试的应用: WB,IP

推荐稀释比: WB: 1/500-1/1000 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R06-0G0

分子量: Calculated MW: 116 kDa; Observed MW: 116 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: -

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TAOK1 in C6, Ramos lysates using TAOK1 antibody.



Western blot analysis of TAOK1 in 3T3, Raw264.7 lysates using TAOK1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
