

TBL1X (4H2) MOUSE MAB

货号: N261005

产品全名: TBL1X (4H2) 小鼠单抗

基因符号 TBL1X; TBL1; F-box-like/WD repeat-containing protein TBL1X; SMAP55; Transducin beta-like protein 1X; Transducin-beta-like protein 1; X-linked

UNIPROT ID: O60907

背景: F-box-like protein involved in the recruitment of the ubiquitin/19S proteasome complex to nuclear receptor-regulated transcription units.

抗原: Purified recombinant human TBL1 protein fragments expressed in E.coli.

经过测试的应用: WB, ICC/IF

推荐稀释比: WB: 1/500-1/1000 IF: 1/50-1/200

种属反应性: Mouse

克隆性: Mouse Monoclonal

克隆编号: 4H2-D5-E9

分子量: Calculated MW: 62 kDa; Observed MW: 58 kDa

亚型: IgG2b


纯化: Affinity Purified

种属反应性: Human


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Neuroscience

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of TBLIX (4H2) in HeLa using TBLIX antibody.



Western blot analysis of TBLIX in SW480 lysates using TBLIX antibody.