

TCEAL6 RABBIT PAB

货号: N225200

产品全名: TCEAL6 兔多抗

基因符号 TCEAL6; Transcription elongation factor A protein-like 6; TCEA-like protein 6; Transcription elongation factor S-II protein-like 6

UNIPROT ID: Q6IPX3

背景: May be involved in transcriptional regulation.

抗原: Synthesized peptide derived from TCEAL6 . at AA range: 1-80

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 22 kDa; Observed MW: 32 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Epigenetics and Nuclear Signaling

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TCEAL6 in CoLo lysates using TCEAL6 antibody.