

## TENSIN 2 RABBIT PAB

货号: N225197

产品全名: Tensin 2 兔多抗

基因符号 TENC1; KIAA1075; TNS2; Tensin-like C1 domain-containing phosphatase; C1 domain-containing phosphatase and tensin homolog; C1-TEN; Tensin-2

**UNIPROT ID:** Q63HR2

背景: Tyrosine-protein phosphatase which regulates cell motility proliferation and muscle-response to insulin.

抗原: The antiserum was produced against synthesized peptide derived from human Tensin-2. AA range:439-488

经过测试的应用: WB,ELISA

推荐稀释比: WB: 1/500-1/1000 ELISA: 1/10000

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 153 kDa; Observed MW: 125 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human,Mouse

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Signal Transduction

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Tensin 2 in A549 3T3 lysates using Tensin 2 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---