

TFG RABBIT MAB

货号: N263028

产品全名: TFG 兔单克隆抗体

基因符号 TF6; HMSNP; SPG57; TRKT3

UNIPROT ID: Q92734

背景: Plays a role in the normal dynamic function of the endoplasmic reticulum (ER) and its associated microtubules.

抗原: A synthetic peptide of human TRK fused gene

经过测试的应用: WB,IHC-P,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Monoclonal

克隆编号: R09-5B1

分子量: Calculated MW: 43 kDa; Observed MW: 55 kDa

亚型: IgG


纯化: Affinity Purified

种属反应性: Human, Mouse and Rat


成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TRK fused gene antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of TRK fused gene in K562, rat Brain, C6, 3T3, Hela lysates using TFG antibody.