

TFG RABBIT PAB

货号: N225732

产品全名: TFG 兔多抗

基因符号 TF6; HMSNP; SPG57; TRKT3

UNIPROT ID: Q92734

背景: Plays a role in the normal dynamic function of the endoplasmic reticulum (ER) and its associated microtubules.

抗原: A synthesized peptide derived from human TRK fused gene

经过测试的应用: WB,IHC-P,ICC/IF,IP

推荐稀释比: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 43 kDa; Observed MW: 55 kDa

亚型: IgG


纯化: Affinity Chromatography

种属反应性: Human, Mouse and Rat

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Immunology

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of TRK fused gene in HeLa lysates using TFG antibody.