

## TGF BETA RECEPTOR I RABBIT PAB

货号: N225337

产品全名: TGF beta Receptor I 兔多抗

基因符号 TGFBR1; ALK5; SKR4; TGF-beta receptor type-1; TGFR-1; Activin A receptor type II-like protein kinase of 53kD; Activin receptor-like kinase 5; ALK-5; ALK5; Serine/threonine-protein kinase receptor R4; SKR4; TGF-beta type I receptor; Transfor

**UNIPROT ID:** P36897

背景: The protein encoded by this gene forms a heteromeric complex with type II TGF-beta receptors when bound to TGF-beta, transducing the TGF-beta signal from the cell surface to the cytoplasm. The encoded protein is a serine/threonine protein kinase.

抗原: Synthetic peptide of human TGFBR1

经过测试的应用: WB,ICC/IF

推荐稀释比: WB: 1/500-1/1000 ICC: 1/100-1/200

种属反应性: Rabbit

克隆性: Rabbit Polyclonal

分子量: Calculated MW: 56 kDa; Observed MW: 56 kDa

亚型: IgG

纯化: Affinity Purified

种属反应性: Human,Mouse


成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

研究领域: Cardiovascular

储存和运输: Store at -20°C. Avoid repeated freezing and thawing



Immunofluorescence analysis of TGF beta Receptor I in Hela using TGF beta Receptor I antibody.



Western blot analysis of TGF beta Receptor I in COLO, HepG2 lysates using TGF beta Receptor I antibody.